



Changsha Sinocare Inc.  
Address: No.265, Guyuan Road, Hi-tech Zone, Changsha, Hunan Province 410205, the People's Republic of China  
Website: [www.sinocare.com](http://www.sinocare.com)  
Email: [intl\\_mkt@sinocare.com](mailto:intl_mkt@sinocare.com)

# Sinocare

Manage Your Blood Glucose with

## iCan i6 CGM



15-day Continuous


Glucose Monitoring



## CGM: a better tool to monitor your glucose

📌 **CGM (Continuous Glucose Monitoring System)** is a system that **continuously** measures glucose levels throughout the **day and night**.

It typically involves a **small sensor** that is **inserted under the skin**, which measures the **glucose in the interstitial fluid** (the fluid between cells).




## Key differences between CGM and BGM

| Feature          | CGM                                       | BGM                                |
|------------------|-------------------------------------------|------------------------------------|
| Monitoring Type  | <b>Continuous, real-time tracking</b> ✓   | Single, discrete tests             |
| Detection Object | <b>Interstitial fluid</b> ✓               | Blood sample (via fingerstick)     |
| Data Frequency   | <b>Every few minutes</b> ✓                | One test at a time                 |
| Alerts           | <b>Yes, for high/low glucose levels</b> ✓ | No alerts; manual testing needed   |
| Comfort          | <b>Small sensor under the skin</b> ✓      | Fingerstick required for each test |
| Duration         | <b>Sensor lasts days to weeks</b> ✓       | Test only at the moment of testing |



## Who can use CGM



**8.71%**

MARD Value

**15** -day

Wear Life 

**Anti-interference** (Acetaminophen and Vitamin C)

Waterproof

**IP28**

2 m 2 h


**iCan**  
CGM System

Use for children over 2 years old



**3** min  
Frequency

**AGP Report 5.0**  
To better help with diabetic complication control



Digital Management



iCan APP  
Patient APP



iCan Reach  
Care Partner APP



iCan Review  
Doctor Software

**Diabetes Management System**

**The world's first 15-day wear life CGM!**

## How to wear Sinocare iCan CGM



Set Up with  
<iCan APP>



Scan the SN Code




Choose the Wearing Site



Avoid areas near waistbands or with scars, tattoos, irritation, and bones.

4

Key Steps



**Disinfect**

Clean the insertion site with **alcohol pad** and allow it to **dry** completely before insertion




**Twist-off**

**Check** the sensor pack before use. Do not use if the label is opened or damaged.



**Insert**


Choose the insertion site and **do not move the applicator**. Press the button in the middle



**Apply Overpatch**

If needed, apply an Overpatch without touching the adhesive area.


## Use iCan APP to manage blood glucose



### Glucose Trend Arrow

- stable glucose trend
- ↘ slow decrease in glucose
- ↓ rapid decrease in glucose
- ↗ slow increase in glucose
- ↑ rapid increase in glucose
- ⇓ sharp decrease in glucose
- ⇓ sharp increase in glucose


## 01/ APP Interface



Record your daily diet, exercise, medication, insulin injection or other events in detail.

## 02/ Alert Settings

- Use the hypoglycemia and hyperglycemia alert settings to **intervene promptly and prevent extreme blood glucose levels.**




### Personalized setting



## 03/ History Review

- By clicking “Events” , users are able to check their History Report over a period of time.



## iCan view Glucose Data Reader: check your glucose without smartphones

- Use the iCan view Glucose Data Reader to display iCan CGM glucose data sent from the transmitter.



### Freedom from Smartphones

Real-time glucose alerts and trends on your standalone Reader via Bluetooth connection.



### Remote Care Synergy

With iCan view WiFi-powered network, families and doctors stay in sync for smarter care.



### Built to Last

IP54 water-resistant, 3-year lifespan.



### Pocket-Sized Power

Fits in pockets, lasts all week.




### Age-Inclusive Design

From toddlers beyond 2 years old to grandparents.



## Share your blood glucosoe data to care partner


- Use the iCan Reach APP to enable remote data sharing and safeguard your blood glucose health.



▲ Invite Care Partners




▲ The Corresponding CGM User Data Viewing




▲ Home Page

## Share your blood glucosoe data to your HCPs


- Use the iCan Review software to enable remote data sharing and safeguard your blood glucose health.



▲ Enter **My Practices** access in CGM application




▲ Input **Personal ID** of HCPs



▲ Data sharing succeeds

## Understand the AGP report

- You can check multi-day glucose data in a 24-hour format AGP report.



## Three steps to assess and improve dynamic blood glucose profiles

- Assess** the dynamic glucose profile



### Monitor Time in Range (TIR)

The longer the time in the target range, the better.



### Assess Blood Glucose Fluctuations

A flatter median curve with a narrower shaded area is preferable.



### Evaluate Hypoglycemia Risk

The shorter the time below range (TBR), the better.

- Improve** the dynamic glucose profile

01

### Reduce Hypoglycemia Risk

Use real-time monitoring to identify periods and causes of hypoglycemia, and consult your doctor to adjust the treatment plan.

02

### Improve Blood Glucose Variability

- ◆ Maintain a regular diet and exercise
- ◆ avoid alcohol, manage stress
- ◆ prevent infections
- ◆ adhere to medication guidelines.

03

### Increase Time in Range

Utilize real-time alerts to intervene promptly, preventing extremes and increasing overall time in the target range.